
DIAION™ Technical Guide

Ion Exchange Resin for Amino Acid Purification

【Amino Acids】

Amino acids are organic compounds made from an amine group and a carboxylic acid group. Amino acid is the smallest unit and component of proteins. Proteins form tissues and organs in human bodies. In addition to acting as a protein component, amino acids also function in biosynthesis inside the body themselves.

Twenty types of amino acids are incorporated into proteins inside the human body. These are divided into essential amino acids, which cannot be synthesized inside the body, and non-essential amino acids. Essential amino acids must be taken through food.

Essential amino acids

L-Isoleucine, L-Leucine, L-Valine, L-Histidine, L-Lysine, L-Methionine, L-Phenylalanine, L-Threonine and L-Tryptophan

Non-Essential amino acids:

L-Asparagine, L-Aspartic Acid, L-Alanine, L-Arginine, L-Cysteine, L-Glutamine, L-Glutamic Acid, Glycine, L-Proline, L-Serine and L-Tyrosine.

Amino acids are industrially prepared by fermentation, chemical synthesis or extraction. Fermentation technology has played crucial roles in this field, and currently the fermented amino acids represent chief products of biotechnology in both volume and value.

The selection of purification methods depends on the nature of fermentation broth.

In case of simple composition broth such as Arginine or Glutamic acid, crystallization is the main process.


If the fermentation broth contains a lot of biomass such as L-Lysine, ion exchange resin method is widely used.

【Lysine Purification】

L-lysine is an essential amino acid for the growth most of animal species and the number one limiting amino acid for poultry. After production and biomass removal by filtration and centrifugation, the essential next step is the lysine purification and recovery. The ion exchange process is one of the most commonly used purification methods.

This process comprises an adsorption step in which a fermentation broth of an amino acid is contacted with a strongly acidic cation exchange resin layer to adsorb the amino acid onto the resin layer. During the adsorption step, suspended substances in the fermentation broth are deposited in the resin so that large amounts of washing water flowing in the reverse direction are needed to remove the deposit after completion of adsorption.

This step is followed by an elution step in which the amino acid is eluted with an eluent such as ammonia water. Eluted L-Lysine is then concentrated and converted into salt form with the addition of acids such as HCL followed by crystallization.

Various processes have been investigated to optimize the production capacity and minimize the water consumption. Most widely used process up to now is a continuous process, in which water and chemical consumption is decreased to one third to half of the conventional process.

【Resin selection for L- Lysine purification】

A wide variety of gel type strong acid cation exchange resins work for most L-lysine recovery operations . A number of issues need to be considered in choosing a resin for the long haul.

1) Type of equipment

The type of equipment influences resin choice.

In case of Large Fixed Bed Systems :Each column is loaded and regenerated over fairly long processing cycles (12 – 24 hours) and a standard gel cation exchange resin such as the Diaion SK1B ,UBK08 can be used in most cases.

In case of a continuous system: resin is used with fairly short processing cycles (2 – 4 hours) and choosing a more osmotically stable and kinetically efficient resin is preferred.

2) Type of broth matrix

Another key consideration is the broth matrix that is being pumped through the resin bed with the L-lysine.

In case of fairly clean broth: Having a settled or filtered broth can help extend resin bed life due to less organic and particulate fouling. Resins with high capacity such as DIAION SKL10 can yield high operation through-put per recovery needs.

In case of broth with an abundance of small peptides: If the broth holds an abundance of small peptides, then choosing a more kinetically open resin, such as DIAION UBK04, offers less potential for fouling over time. Mitsubishi offers a wide range of cation exchange resins for your L-lysine recovery needs. The attached table is a good starting point.

Product	Description	Type	Structure	Matrix	Total Exchange Capacity
DIAION™ SK 1 B	Osmotically strong, gel type standard cation exchange resin	Strong acid cation exchange resin	Gel	Sulfonated Styrene/DVB copolymer 8% crosslinked	Min, 2mol/L
DIAION™ SKL10	Osmotically strong, high capacity, uniform particle size gel type cation exchange Resin	Strong acid cation exchange resin	Gel	Sulfonated Styrene/DVB copolymer 10% crosslinked	Min. 2.2mol/L
DIAION™ UBK08	Osmotically strong, uniformed particle size gel type cation exchange resin	Strong acid cation exchange resin	Gel	Sulfonated Styrene/DVB copolymer 8% crosslinked	Min. 2 mol/L
DIAION™ UBK04	Kinetically fast, uniformed particle size gel type cation exchange resin	Strong acid cation exchange resin	Gel	Sulfonated Styrene/DVB copolymer 4% crosslinked	Min. 1. 2 mol/L

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DIAION™ and SEPABEADS™ is a registered trademark of Mitsubishi Chemical Corporation.
RELITE™ is a registered trademark of Resindion S.R.L.
Mitsubishi Chemical Aqua Solutions Co., Ltd. (Tokyo)
Tai-Young Chemical Co, Ltd. (Kaohsiung)
Mitsubishi Chemical (China) Co., Ltd (Shanghai)
Mitsubishi Chemical India Pvt. Ltd. (Haryana)