


# DIAION™ SEPABEADS™ MCIGEL™ Technical Guide

## Polymeric Separation Media for Pharmaceutical Purification

Separation of 6-amino-penicilinic acid derivatives

Separation of cyclic peptide analogs

Mitsubishi Chemical's polymeric separation media have been used in pharmaceutical application from late 1960's and have become one of the standard resins in purification not only in pharmaceuticals but also juice, sugar and food segment.

Mitsubishi Chemical offers wide variety of polymeric separation media under the trade name "DIAION", "SEPABEADS" grades for industrial application and "MCIGEL CHP", "MCIGEL CMG" grades for analytical and preparative applications.

These media have been adapted in the separation/purification process of antibiotics, oligonucleotide, peptides, herb and Chinese medicines.

# Separation and Purification of Pharmaceuticals

Peptides and nucleic acids used as drug substances for pharmaceuticals are produced by chemical synthesis (solid-phase or solution) or by recombinant procedures. In both cases, reversed-phase chromatography (RPC) is a common method for separation and purification of target molecules.

Mitsubishi Chemical's polymeric separation media such as DIAION HP20, SEPABEADS SP20SS and MCIGEL CHP20 series are commonly used as stationary phase of RPC process for industrial-scale pharmaceutical purification.

## Principle of Reversed-Phase Chromatography (RPC)

In RPC, substances in mobile phase tend to be retained on stationary phase in chromatography column via hydrophobic interactions. The strength of the interaction generally depends on hydrophobicity of the stationary phase and substances and relationship between pore size of the stationary phase and molecular size of the substances. Substances with strong interaction have a longer retention time in the column and the difference in retention time leads to separation of each substance.

Polymeric separation media have covalently bound hydrophobic moieties such as aromatics in their chemical structure and are suitable as stationary phase for RPC.



Illustration of principle of RPC

## General procedure and equipment of RPC

The general procedure of RPC is as follows.

The solution that contains mixture of target molecules and impurities is loaded to the column packed with separation media (stationary phase), followed by elution with eluting solvent (mobile phase). The mobile phase should be able to dissolve the target molecules. For industrial-scale pharmaceutical purification process, aqueous solution of acetonitrile or alcohol is commonly used as mobile phase.

The schematic equipment used in RPC is shown in the figure. Component molecules in the feed solution travel through the stationary phase at different speeds depending on their individual natures and elute from the column separately. Elution of the molecules from the column is detected by UV detector etc. and purified target molecule is selectively fractionated.

Chromatographic separation equipment

# Polymeric Separation Media for Pharmaceutical Purification

Mitsubishi's polymeric separation media are the spherical porous resins based on crosslinked polystyrene, crosslinked modified polystyrene or crosslinked methacrylate.

The advantages of the products can be summarized as follows;




- High chemical stability over all pH range
- High physical strength
- Excellent cycle stability and long life
- Wide variety of degree in hydrophobicity based on polymer chemistry
- Wide variety of particle size from preparative to industrial application

## How to select product grade

Mitsubishi's polymeric separation media can be roughly divided into three types based on the composition of base polymer resin. Depending on the nature of target molecule, the right product can be generally selected as follows.

- DIAION HP20, SEPABEADS SP20SS and MCIGEL CHP20 series are the standard styrene-DVB resins with moderate hydrophobicity, so these can be selected as a first choice for new process.
- SEPABEADS SP207 and other brominated-styrenic resins have higher hydrophobicity which is suitable for purification of high polarity molecules. Also, these grades have high specific density and are easy to handle.
- DIAION HP2MGL and other methacrylic polymer resin have less hydrophobicity than styrenic resins and are suitable for purification of highly hydrophobic molecules. These can also be used for purification of oil soluble molecules using normal phase chromatography.


Also, each type of product has a variety of particle size grades and can be selected according to the balance between separation performance and pressure restriction.

Polymer structure	Resin grades	Particle size (μm)	Application (extraction, separation, purification)
<b>Polystyrene</b> (hydrophobic)  	HP20,HP21 SP825L,SP850 SP70, SP700	>250	Antibiotics (Cephalosporin, Vancomycin, Carbapenem and others.) Anti-cancer drug, immunosuppressive drug
	HP20SS SP20SS	120 70	Antibiotics (Vancomycin, Daptomycin) Anti-cancer drug, immunosuppressive drug
	CHP20 series	20,30,50	Peptide, Oligo-nucleotide
<b>Modified Polystyrene</b> (Highly hydrophobic)  	SP207	>250	Antibiotics (Cephalosporin and others) Low hydrophobic drug
	SP207SS	120	Low hydrophobic drug
	CHP07 (to be launched)	30	Peptide, Oligo-nucleotide
<b>Methacrylate</b> (Low hydrophobic)  	HP2MGL	>250	Cyclic Peptide (anti-cancer drug) Highly hydrophobic drug
	HP2MGSS (to be launched)	120	Highly hydrophobic drug Oil soluble drug (Normal phase chromatography)
	CMG20	30	Peptide, Oligo-nucleotide

# Product Selection Guide


## • Selection by nature of target molecule

The type of separation media can be selected based on the nature of target molecule to be separated. In general, highly hydrophobic media are suitable for separation of low hydrophobic (high polarity) molecules, and low hydrophobic media are suitable for highly hydrophobic molecules.


## Comparison of retention time of cyclic peptide using polymeric media with different hydrophobicity

Under the same analytical conditions, the separation media with higher hydrophobicity shows longer retention time.

CF: C<sub>35</sub>H<sub>54</sub>N<sub>6</sub>O<sub>6</sub>  
MW: 654.8530

Column: 4.6 mm I.D. x 250 mm  
Temperature: 35°C  
Flow rate: 1 mL/min  
Solvent: 65-95% CH<sub>3</sub>CN in H<sub>2</sub>O, 20 min linear gradient



## • Selection by particle size

The particle size of separation media can be selected based on the required separation performance and pressure characteristics. In general, small particles are suitable for precision purification due to high separation performance, and large particles are suitable for large-scale refining due to low pressure.

	Small particle (20 - 50 $\mu\text{m}$ )	Middle particle (70 - 120 $\mu\text{m}$ )	Large particle (over 250 $\mu\text{m}$ )
<b>Polystyrene</b> (hydrophobic)	CHP20 (20, 30, 50 $\mu\text{m}$ )	SP20SS (70 $\mu\text{m}$ ) HP20SS (120 $\mu\text{m}$ )	HP20, HP21. SP825L, SP850 SP70, SP700
<b>Modified Polystyrene</b> (Highly hydrophobic)	CHP07 (30 $\mu\text{m}$ ) (to be launched)	SP207SS (70 $\mu\text{m}$ )	SP207
<b>Methacrylate</b> (Low hydrophobic)	CMG20 (30 $\mu\text{m}$ )	HP2MGSS (120 $\mu\text{m}$ ) (to be launched)	HP2MGL

## Separation of penicillanic derivatives using polystyrenic media with different particle size

The chromatographic separation data below shows the effect of particle size on the separation of three penicillanic derivatives using polymeric media with same chemistry. The smaller particle shows sharper peaks and better separation performance.


Column size: 250 x 10mm I.D., Eluent: MeOH/50mM phosphate (pH8.0) = 60/40  
 Flow rate: 2.18ml/min, Detection: UV 254nm  
 Samples: 1: 6-aminopenicillanic acid (1g/L)  
           2: penicillin G (1g/L)  
           3: penicillin V (1g/L).  
 Injection: 100 $\mu\text{L}$

# Examples of Separation using Polymeric Separation Media


## Separation of cyclic peptide analogs

CHP20 can separate cyclic peptides with six similar chemical structures as shown below. The separation is carried out using acetonitrile/water gradient as eluent.


## Sample loadability of polymeric separation media


Compared to ODS-silica media, CHP20 maintains good separation behavior even at high loadings in cyclosporine analogs separation.

**Low loading condition** (Sample injection: 0.01 mg each of CsA, CsB, CsC, and CsD)

**High loading condition** (Sample injection: Total 0.1, 0.2, 0.5 and 1.0 mg (content ratio: CsA 90%, CsB 3.3%, CsC 3.3%, CsD 3.3%))

Column: 4.6 mm x 250 mm, Temp: 35°C, Flow rate: 1 mL/min, Detection: 214 nm  
Solvent: (CHP20) 69% CH<sub>3</sub>CN in H<sub>2</sub>O, 30 min isocratic elution, (ODS) 75% CH<sub>3</sub>CN in H<sub>2</sub>O, 30 min isocratic elution

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DIAION™, SEPABEADS™ and MCI GEL™ is a registered trademark of Mitsubishi Chemical Corporation.
RELITE™ is a registered trademark of Resindion S.R.L.
Mitsubishi Chemical Aqua Solutions Co., Ltd. (Tokyo)
Tai-Young Chemical Co, Ltd. (Kaohsiung)
Mitsubishi Chemical (China) Co., Ltd (Shanghai)
Mitsubishi Chemical India Pvt. Ltd. (Haryana)